


**NEW**



**TOPIC:**

**SAMPLE TYPE**

**URINE**  
**(undiluted)**



# Additional sample type: undiluted urine



## NOW AVAILABLE



no sample preparation




no additional reagents





# Measured Parameters



		Range	CV in % <sup>*</sup>
<b>Na<sup>+</sup></b>	<b>Sodium</b>	5.0 - 300 mmol/l	0.50
<b>K<sup>+</sup></b>	<b>Potassium</b>	5.0 - 120 mmol/l	0.73
<b>Cl<sup>-</sup></b>	<b>Chloride</b>	5.0 - 300 mmol/l	0.62

<sup>\*</sup> **Repeatability** (within-run precision)



11:39

# Cartridge information > Details

Cartridge state ✓ valid

		Status
Time since last fluidic action	3 h	✓
Urine samples	99 / 150	✓

← Back



How many  
urine  
samples?

150 samples / cartridge

**WE CHALLENGE  
THE STATUS QUO!**

